



PEAKS Online: Automated Quality Control (QC) Tool

Yiming Xiao, PhD, Application Scientist
Kyle Hoffman, PhD, Applications Manager

Bioinformatics Solutions Inc., Waterloo, Canada


Aim:

To demonstrate the utility of the Automated Quality Control (QC) tool PEAKS Online 11.

Introduction:

LC-MS based protein/peptide quantification has drawn attention for physiopathology/pharmaceutical studies. For such research, especially clinical studies, obtaining a large number of samples is necessary to provide sufficient statistical power, i.e., over 50 samples per group are often necessary due to the high interindividual variation [1]. For cohorts with large sample sizes, performing a QC step is essential for ensuring data quality and validity of the results.

PEAKS Online has a automated QC tool that provides sophisticated QC analysis on top of protein/peptide identification and quantification results. All the statistical analyses are user-friendly and can be specified based on the user's requirements.



Case study: benchmarking data set analysis with PEAK Online:

A published data set [2] was used as an example to demonstrate an application of LFO and QC analysis in PEAKS Online.

Study aims and background:

Label-free quantification using DDA and DIA (SWATH) can be compared to determine which acquisition method leads to more reproducible and accurate results. DDA is the more traditional way but also suffers from high missing values and under-sampling. DIA is more robust and can resolve DDA-related biases, however, DIA also has issues such as high false positive rate in complex MS2 spectra. In this study, the same samples were analyzed by high-resolution DDA or DIA (SWATH) and compared. When performing MS1 ion current extraction, high resolution DDA (HS-DDA) quantification was comparable to DIA in accuracy, precision, and better for lower abundance proteins. In this application note, the DDA data is used for presenting automated QC tool in PEAKS Online, to highlight attributes that fail in the analysis.





b Match Between Runs

Mass Error Tolerance: Tolerance Unit: Retention Time Shift Tolerance(min): Auto Detect

Feature Intensity \geq

RT Range: \leq RT \leq Base Sample:

Peptide feature Avg. Area $\geq 2e+5$, Quality ≥ 20 , $2 \leq$ Charge ≤ 5 , Peptide ID Count ≥ 1 per a group, and detected in at least 1 samples per a group.

Protein Significance Method: ANOVA, Modified Form Exclusion, Remove Outlier, Use Top 3 peptide, Significance ≥ 0 , $1 \leq$ Fold change ≤ 64 , has at least 1 used peptide

Normalization Method: Use TIC

Fig 1. Protein/peptide quantification in PEAKS Online. (a) PEAKS Q workflow integrated with QC function. (b) Lfq parameters used.

Methods:

Five groups of samples were prepared with three different proteomes (human, E. coli and yeast) and run with five technical replicates (n=25). The human protein proportion was 60% across all samples. The portion of E. coli to yeast protein amounts were as follows:

A: 5%/35%, B: 7.5%/32.5%, C: 10%/30%, D: 15%/25%, E: 20%/20%.

DDA data was acquired with MS1 240000 resolution and MS2 15000 resolution. The gradient was 160 minutes long. A detail of LC-MS method can be found in [2].

Results:

MS data (5 runs *5 samples) was analyzed in PEAKS Online as Lfq with PEAKS Q module (Figure 1a). For the quantification part, match-between-run and TIC normalization were applied. The detailed search parameters are shown in Figure 1b.




Fig 2. Lfq result. (a) Number of identified attributes. (b) The volcano plot for differentially expressed proteins.

While using at least two peptides per protein for quantification, 6076 protein groups are quantifiable. The number of features, features with identifications, and protein groups are listed in Figure 2a. The volcano plot shows the differential expression of proteins across all samples (Figure 2b).




Fig 3. (a) Number of features found in each sample. (b) Cumulative precursor count plotted against retention time.

Results cont'd:

QC result of LC-MS data

In QC analysis, the sample average was set as the reference, and the acceptance tolerance was set as 10%. Any attribute that falls outside of the 10% acceptance tolerance is labelled in red (fail) in the QC result views (Table 1).

The QC result shows the number of MS1, MS/MS, features, MS2/MS1 ratio, full peak width, full peak width at half maximum (FWHM), and total base peak chromatogram (BPC) intensity. The data QC result shows that out of all 25 samples (E01-E05 are not shown in Table 1) and 7 different metric categories, only B04 and B05 have a lower BPC intensity compared to the average value (12.7% and 11.6%), causing this attribute to fail and fall outside of the 10% acceptance tolerance.

← 1 - 21 of 21 →

| SAMPLE NAME | # MS1 | # MS/MS | MS2/MS1 RATE | # FEATURES | FULL WIDTH (MIN) | FWHM (SEC) | TOTAL BPC |
|-------------|-------|---------|--------------|------------|------------------|------------|-----------|
| Average | 5773 | 72711 | 12.60 | 230574 | 0.31 | 11.21 | 3.882e+9 |
| A01 | 5641 | 73581 | 13.04 | 239956 | 0.31 | 11.31 | 4.022e+9 |
| A02 | 5882 | 71706 | 12.19 | 232641 | 0.31 | 11.23 | 3.754e+9 |
| A03 | 5874 | 71802 | 12.22 | 229834 | 0.31 | 11.22 | 3.890e+9 |
| A04 | 5825 | 72155 | 12.39 | 229729 | 0.31 | 11.21 | 3.737e+9 |
| A05 | 5871 | 71830 | 12.23 | 221794 | 0.31 | 11.18 | 3.664e+9 |
| B01 | 5684 | 73493 | 12.93 | 237093 | 0.31 | 11.20 | 4.060e+9 |
| B02 | 5792 | 72458 | 12.51 | 232858 | 0.31 | 11.20 | 3.867e+9 |
| B03 | 5799 | 72375 | 12.48 | 232093 | 0.31 | 11.20 | 3.718e+9 |
| B04 | 5778 | 72370 | 12.53 | 224801 | 0.31 | 11.19 | 3.389e+9 |
| B05 | 5807 | 71887 | 12.38 | 228044 | 0.31 | 11.19 | 3.431e+9 |
| C01 | 5680 | 73690 | 12.97 | 236578 | 0.31 | 11.20 | 4.172e+9 |
| C02 | 5787 | 72654 | 12.55 | 232673 | 0.31 | 11.20 | 3.960e+9 |
| C03 | 5834 | 72219 | 12.38 | 231485 | 0.31 | 11.20 | 3.824e+9 |
| C04 | 5834 | 72473 | 12.42 | 225205 | 0.31 | 11.20 | 3.782e+9 |
| C05 | 5863 | 71803 | 12.25 | 227010 | 0.31 | 11.19 | 3.691e+9 |
| D01 | 5661 | 73921 | 13.06 | 235646 | 0.31 | 11.20 | 4.235e+9 |
| D02 | 5819 | 72415 | 12.44 | 231106 | 0.31 | 11.20 | 4.034e+9 |
| D03 | 5783 | 72649 | 12.56 | 228636 | 0.31 | 11.20 | 3.899e+9 |
| D04 | 5724 | 72882 | 12.73 | 229480 | 0.31 | 11.21 | 3.696e+9 |
| D05 | 5775 | 72734 | 12.59 | 227576 | 0.31 | 11.21 | 3.845e+9 |

Table 1. LC-MS Data QC Result




Fig 4. Sample TIC over retention time for samples A01 to A05.

The number of features (highlighted in blue in Table 1) in each sample is also presented in a bar chart with min and max tolerances (Figure 3a). The cumulative precursor count is shown in Figure 3b. A near linearity increment of precursor across retention time suggests the peptides were well separated by the LC gradient.

The Chromatogram of each sample (up to 5) overlay in Figure 4, provides an overview of total ion chromatograms (TIC) from the selected LC-MS runs. This visualization allows the user to assess differences in the chromatograms across sample runs.

| SAMPLE NAME | #MS/MS | # IDENTIFIED FEATURES | # PSM | # IDENTIFIED PEPTIDES | # QUANTIFIED PEPTIDES | # SEQUENCES | # IDENTIFIED PROTEIN GROUPS | # QUANTIFIED PROTEIN GROUPS | # TOP PROTEINS | # ALL PROTEINS | # SCANS / # PEPTIDE RATE | ID RATE |
|----------------|--------|-----------------------|-------|-----------------------|-----------------------|-------------|-----------------------------|-----------------------------|----------------|----------------|--------------------------|---------|
| Average | 72711 | 47774 | 59700 | 42157 | 36387 | 41565 | 7049 | 6064 | 7139 | 7493 | 1.22 | 82.10% |
| A01 | 73581 | 50751 | 62123 | 44946 | 36424 | 44402 | 7088 | 6066 | 7179 | 7564 | 1.18 | 84.43% |
| A02 | 71706 | 46906 | 58171 | 41330 | 36229 | 40688 | 6794 | 6049 | 6881 | 7219 | 1.23 | 81.12% |
| A03 | 71802 | 46686 | 58623 | 41474 | 36386 | 40873 | 6908 | 6067 | 6994 | 7349 | 1.22 | 81.65% |
| A04 | 72155 | 47691 | 59931 | 42354 | 36400 | 41754 | 7067 | 6065 | 7160 | 7520 | 1.20 | 83.06% |
| A05 | 71830 | 45214 | 58250 | 39899 | 36236 | 39415 | 6793 | 6053 | 6885 | 7223 | 1.23 | 81.09% |
| B01 | 73493 | 49832 | 61447 | 44158 | 36438 | 43544 | 7159 | 6072 | 7246 | 7621 | 1.20 | 83.61% |
| B02 | 72458 | 48131 | 59619 | 42513 | 36356 | 41845 | 7089 | 6062 | 7180 | 7530 | 1.22 | 82.28% |
| B03 | 72375 | 47061 | 58603 | 41664 | 36443 | 41066 | 6938 | 6071 | 7016 | 7356 | 1.24 | 80.97% |
| B04 | 72370 | 45055 | 57661 | 39567 | 36275 | 39050 | 6842 | 6057 | 6926 | 7248 | 1.26 | 79.68% |
| B05 | 71887 | 46539 | 58916 | 41184 | 36368 | 40629 | 6924 | 6057 | 7012 | 7378 | 1.22 | 81.96% |
| C01 | 73690 | 49720 | 61302 | 43866 | 36412 | 43217 | 7142 | 6062 | 7232 | 7598 | 1.20 | 83.19% |
| C02 | 72654 | 47421 | 58909 | 41880 | 36418 | 41267 | 7044 | 6065 | 7131 | 7479 | 1.23 | 81.08% |
| C03 | 72219 | 47622 | 59532 | 42228 | 36441 | 41617 | 7101 | 6072 | 7182 | 7537 | 1.21 | 82.43% |
| C04 | 72473 | 47192 | 59254 | 41721 | 36294 | 41172 | 7190 | 6062 | 7275 | 7632 | 1.22 | 81.76% |
| C05 | 71803 | 46176 | 58616 | 40828 | 36413 | 40269 | 6908 | 6068 | 6986 | 7331 | 1.22 | 81.63% |
| D01 | 73921 | 49407 | 60493 | 43440 | 36413 | 42828 | 7205 | 6069 | 7311 | 7688 | 1.22 | 81.83% |
| D02 | 72415 | 48191 | 60110 | 42553 | 36415 | 41941 | 7131 | 6061 | 7218 | 7579 | 1.20 | 83.01% |
| D03 | 72649 | 47558 | 59381 | 42138 | 36420 | 41541 | 7120 | 6072 | 7208 | 7558 | 1.22 | 81.74% |
| D04 | 72882 | 47452 | 59553 | 41690 | 36438 | 41122 | 7109 | 6072 | 7206 | 7561 | 1.22 | 81.71% |
| D05 | 72734 | 47235 | 59792 | 41503 | 36427 | 40917 | 7014 | 6066 | 7099 | 7459 | 1.22 | 82.21% |

Table 2. LFQ QC Result.

QC result of LFQ

Similar to the QC result of LC-MS data, a table is provided for the LFQ QC result (Table 2). The highlighted part (identified and quantified protein groups) is also presented in a bar chart (Figure 5a). The min and max acceptance tolerance is calculated based on the average number of quantified protein groups. The number of proteins with the corresponding covariance values (CV) are shown in Figure 5b. The CV value here reflects the variance in a protein's abundance across all samples.




Fig 5. Number of identified and quantified protein groups and covariance values. (a) Identified and quantified protein groups after match between runs. (b) Number of proteins with corresponding CV values.

The peptide/protein identification reproducibility figures (Venn diagrams) shown in Figure 6 allow the user to quickly determine how many common or unique peptides/proteins are identified between samples. The Pearson correlation charts of samples shown in Figure 7 displays the correlation of peptide or protein intensities between the two selected samples. The user could pick any pair of samples to perform such examinations.




Fig 6. Sample reproducibility between A01 and A02. (a) Number of peptides unique or common in samples A01 and A02. (b) Number of proteins unique or common in samples A01 and A02.




Fig 7. Pearson correlation of A01 and A02. (a) Correlation of peptides in samples A01 and A02. (b) Correlation of proteins in samples A01 and A02.




Fig 8. Box plot of the difference in the retention times of peptides from all samples across the LC gradient, with respect to a specified control sample.

Conclusion:

The new version of PEAKS Online integrates an automated QC tool for protein/peptide identification and quantification. By coupling QC information with data analysis, PEAKS Online can help users efficiently validate sample data in large cohorts and provide information for any potential troubleshooting.

References:

1. Overmyer, K. A. et al. Large-scale multi-omic analysis of COVID-19 severity. *Cell Syst.* 12, 23–40 (2021).
2. Wang, X. et al. Ultra-High-Resolution IonStar Strategy Enhancing Accuracy and Precision of MS1-Based Proteomics and an Extensive Comparison with State-of-the-Art SWATH-MS in Large-Cohort Quantification. *Analytical Chemistry.* 11, 4884-4893 (2021).

Bioinformatics Solutions, Inc.

140 Columbia St, Suite 202
Waterloo, Ontario N2L 3K8
Canada

Tel: (519) 885-8288
Fax: (519) 885-9075

sales@bioinfor.com
www.bioinfor.com



Information, descriptions, and specifications in this publication are subject to change without notice.
Bioinformatics Solutions, Inc. 2023